



Oxysure Go Portable Oxygen Concentrator



1.6kg
Excl. Battery



1-5
Pulse Flow



Upto 10 hrs
(Double Battery)



Weight

1.6 kg

0.5 kg per battery extra

Reliable Battery for Everyday Use

The Oxysure Go is designed to support your daily routine without frequent recharging. Whether you are at home, outdoors, or travelling, the long-lasting battery gives you confidence and freedom throughout the day.

Pulse Setting	Duration (Single Battery)	Duration (Double Battery)
1	5 h 00 min	10 h 00 min
2	3 h 50 min	7 h 40 min
3	3 h 00 min	6 h 00 min
4	2 h 00 min	4 h 00 min
5	1 h 40 min	3 h 20 min

*under ideal conditions – actual backup may vary based on operating conditions and breath rate

Recharge time: 2 hours per battery



Designed for Comfort and Portability

Its compact size fits comfortably in daily life easy to carry, store, and travel with, without feeling bulky or heavy.

Clean, Safe Oxygen You Can Trust

Oxysure Go delivers clean, reliable oxygen using an **oil-free compressor**, ensuring **no risk of contamination**. It is quiet, dependable, and built for long-term daily use.



Button Bump Design

Button bump design makes controls easy to touch and operate.



Push-Pull Battery Design

One-button battery removal, simple to use and elderly-friendly.




Mushroom Head Oxygen Outlet

Prevents oxidation related cracking and securely holds oxygen tubing in place.

Intuitive Panel and Display

Easy-to-read screen and intuitive buttons ensure confident, stress-free daily use.



Move freely, with reliable oxygen by your side.

Oxysure™ Go is designed for active lifestyles, delivering dependable oxygen wherever life takes you. Quality-tested and easy to carry, it supports your freedom to move.

Oxysure Go Specifications

Weight	1.6kg (Machine), 0.5kg (Battery)
Size	Length: 183mm, Width: 86mm, Height: 199mm
Oxygen Flow	Pulse flow delivery system, 1-5 pulse settings Setting 1: 210mL (Equivalent to 0.6LPM) Setting 2: 420mL (Equivalent to 1.2LPM) Setting 3: 640mL (Equivalent to 1.9LPM) Setting 4: 840mL (Equivalent to 2.5LPM) Setting 5: 1000mL (Equivalent to 3LPM) <i>Equivalent flow output may vary between patients and is calculated based on an I:E ratio of 1:2. Please consult your physician and monitor SpO₂ after setting the flow.</i>
Sound Level	<49dB, Quite Operation
Purity	93% + 3% at all settings
Airline Use	Meets FAA acceptance criteria for carriage and use on board aircraft
FDA	Approved
Warranty	3 years (Main unit), 1 year (Battery & sieve beds)
Operating Altitude	Upto 13,123 feet (4000m)
AC Power	100-240V (50-60Hz), Auto-Sensing Worldwide Power
Environment Range	-10 to 40°C, Humidity: <80%, Atmospheric Pressure: 86-106Kpa
Battery Duration	Up to 10 hours (with Double Battery)
Recharge Time	2 hours per battery



Disclaimer: This portable oxygen concentrator provides supplemental oxygen only and is not a life-support device. Use strictly as prescribed by a medical professional. Oxygen delivery and battery performance may vary based on usage conditions. Not intended for emergency or critical care use. Airline acceptance is subject to individual airline policies; users are advised to confirm requirements with the airline prior to travel. Specifications subject to change without notice.