

Oxymed P2

Portable Oxygen Concentrator















Effortless Portability




Enjoy the Freedom of Portable Oxygen Therapy, enabling you to travel, shop, and engage in activities with ease and peace of mind.





Easy Control system with buttons and intuitive design

P2 offers effortless navigation with user-friendly buttons, putting ultimate relaxation control right at your fingertips






Approved for In-flight use

Drive with Oxygen

Keep your life-saving oxygen concentrator charged while driving, ensuring a constant supply of breathable air



Impressive Battery Backup up to 10 hours

Never worry about running out of power with our exceptional battery backup, giving you peace of mind for up to 10 hours of usage.

Pulse Setting	Duration (Single Battery)	Duration (Double Battery)	
1	5 h 00 min	10 h 00 min	
2	3 h 50 min	7 h 40 min	
3	3 h 00 min	6 h 00 min	
4	2 h 00 min	4 h 00 min	
5	1 h 40 min	3 h 20 min	

Recharge Time: 4 hours per battery

Uninterrupted PowerUse one, Charge the other

External Fast Charger*

Conveniently charge a backup battery while the device is in use.

Charging Time: 2.5 hours *Additional purchase.



Oxymed P2 Portable Oxygen Concentrator

Specifications

Weight (with single battery)	1.98 kg	Power Supply	AC Input: 100 to 240VAC 50 to 60 Hz DC Input: 11-16 VDC
Size	8.70"L x 6.30"W x 3.35"H (includes battery)		DC Output: 19V 6.3A
Noise Level	49 dB @ Setting 2	Warranty & Quality Standards	
Use	Designed for Continuous Use	Quality Standards	FDA, CE, FAA Approved
Flow settings & Pulse Volume		In-Flight Use	Approved
Oxygen Flow Settings (Pulse)	Pulse Dose Delivery System Upto 5	Warranty	3 years (Main unit) 1 year (Battery and sieve beds)
Maximum Minute Volume Output	Setting 1 - 210 mL / min Setting 2 - 420 mL / min Setting 3 - 630 mL / min Setting 4 - 840 mL / min Setting 5 - 1000 mL / min	What's in the Box	Oxymed P2 unit AC Adapter Car Charger Nasal Cannula Single Battery Sling bag Screwdriver 1 set of filters User manual Double battery is an additional purchase
Oxygen Concentration	90% - 3%/+6% at all settings		
Warm up Time	2 minutes		
Display Screen	2.8 inch	-	
Battery & Power			,,
Battery Backup	Single Battery (Included) Upto 5 hours	-	
	Double Battery (Additional Purchase) Upto 10 hours		





Registered Office:

A brand by Medequip Healthcare Solutions LLP Plot No: A-12, A-13, Apparel Park Phase-1, Doddaballapur Industrial Area Bengaluru - 561203, Karnataka

Branch Offices:

Delhi, Telangana, Odisha, Gujarat, Maharashtra, Uttar Pradesh, Rajasthan, Madhya Pradesh Punjab, Tamil Nadu, West Bengal, Kerala, Jharkhand, Chhattisgarh and Bihar





