




Oxymed i Series B5





A sleep therapy solution for Non complaint Obstructive Sleep Apnea patients (OSA) with OHS (Obesity Hypoventilation Syndrome)

Oxymed i series B5 Auto BPAP is designed to help and manage your challenging sever OSA patients with OHS and also helps to treat patient with overlap syndrome. A new i series B5 device offers a new beginning in a great sleep.

Effective Therapy

- Advanced Algorithm for Central apnea detection
- AST technology for better leak management
- Data management
7 days of detailed data and 30 days of therapy data,
180 days of summary data, with 16 GB pen drive storage
- Quiet Therapy

Easy to use

- Compact and Smart
- Auto On/Off
- One button control
- preheat
- Intuitive user interface with an easy-to-read angled screen

Greater Comfort

- COMF (pressure release technology)
- Ramp OFF/0-45 minutes (5 minutes interval)
- Fully Integrated Humidifier
- High efficiency blower





Technologies working together.

AST (Automatic Synchronization Technology)

This feature is designed to adjust the pressure based on mask leak, to make sure the set pressure is delivered to patient,

Comfort features

One button control

Easy to operate the device, and we can enable the lock settings - if in case someone doesn't want to disturb the settings.

Preheat

It helps user to get the optimal temperature from the 1st minute of device usage, device can rise the temperature in 5 minutes from 23°C to 45°C.




Integrated Humidifier

Oxymed I series B5 comes with integrated humidifier, so the patient can enjoy the humidification. Heated humidification may help to reduce the common side effects of CPAP/BIPAP therapy like Nasal Congestion, dry nose during winter etc.

COMF (Pressure Relief Technology)

Advanced and expiratory pressure relief function/COMF. It has settings of 0 & 1-3 levels, helps patient to release or exhale CO₂ freely. Basically, COMF will help to reduce the pressure. It delivers comfortable breathing experience to the user.



I series BIPAP (New edition)

Sleep Apnoea Therapy Devices

Product Specifications

Model	B5
Mode	CPAP, AUTO CPAP, BPAP, AUTO BPAP
Pressure Range	4-30 cm H ₂ O
Pressure Accuracy	+/-2cm H ₂ O
Ramp Time	0 min-45min(5-minute increments) , Ramp pressure: 4cm H ₂ O~ (EPAP-0.5) cmH ₂ O
COMF Pressure Relief	1-3 levels
Humidification Level	1-5 levels (95 to 167°F/35 to 75°C)
BPM	1-30
I:E (S,S/T)	1:3
Parameter Locking	Yes
Auto Drying	Yes
rise time	1-6 levels 1-4 levels
Data storage capacity	16GB
Sleep Data Memory	180 days
Noise	32 dB (A) Under the condition that the pressure is 10cmH ₂ O in CPAP mode)



REGISTERED OFFICE :-

MEDEQUIP HEALTHCARE SOLUTIONS LLP,
PLOT No: A-12, A-13, APPAREL PARK,
PHASE-1, DODDABALLAPUR INDUSTRIAL AREA,
BANGALORE-RURAL- 561203, KARNATAKA.

Branch Offices:

Delhi, Telangana, Odisha, Gujarat, Maharashtra
Uttar Pradesh, Rajasthan, Madhya Pradesh
Punjab, Tamilnadu, West Bengal, Kerala
Jharkhand, Chhattisgarh and Bihar

To stay uptodate about
our products, follow us on :



080-28362444, 080-28362555



www.medequip.co.in